
Product Name: CD81 (1803) Rabbit Monoclonal Antibody**Catalog #: AMRe08464**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	IHC,ICC/IF,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	IHC 1:200-1:1000,ICC/IF 1:100-1:200,FC 1:200-1:500
Molecular Weight	26kDa

Antigen Information


Gene Name	CD81
Alternative Names	
Gene ID	975.0
SwissProt ID	P60033
Immunogen	Recombinant protein of human CD81

Background

Research Area

Tags & Cell Markers

Image Data



Western blot detection of CD81 in C2C12, 3T3, MCF-7 cell lysates using CD81 antibody(1:1000 diluted).