

**Product Name: Cdc6 (14X12) Rabbit Monoclonal Antibody****Catalog #: AMRe08528**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.5mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Purification</b>	Affinity purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:100-1:200
<b>Molecular Weight</b>	63kDa

**Antigen Information**

<b>Gene Name</b>	CDC6
<b>Alternative Names</b>	CDC6; CDC18L; CDC6-related protein; Cdc18-related protein; HsCDC6; HsCDC18; p62(cdc6);
<b>Gene ID</b>	990.0
<b>SwissProt ID</b>	Q99741
<b>Immunogen</b>	A synthetic peptide of human Cdc6

**Background**


Involved in the initiation of DNA replication. Also participates in checkpoint controls that ensure DNA replication is completed

before mitosis is initiated. Cdc6 has recently been shown to play an important role in the intra-S-phase p21 Waf1/Cip1-dependent DNA damage response. Both cdc6 and CDT1 are degraded by the ubiquitin proteasome pathway in response to DNA damage associated with re-replication. Involved in the initiation of DNA replication. Also participates in checkpoint controls that ensure DNA replication is completed before mitosis is initiated.

## Research Area

Epigenetics and Nuclear Signaling

## Image Data



Western blot analysis of Cdc6 expression in HeLa cell lysate treated with hydroxyurea.