

Product Name: Claudin 11 (15R12) Rabbit Monoclonal Antibody
Catalog #: AMRe08891



Summary

Production Name	Claudin 11 (15R12) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	CLDN11
Alternative Names	Claudin 11; CLD11_HUMAN; Cldn11; OSP; OTM;
Gene ID	5010.0
SwissProt ID	O75508.

Application

Dilution Ratio	WB 1:1000-1:5000, FCM 1:100-1:500
Molecular Weight	22kDa

Product Name: Claudin 11 (15R12) Rabbit Monoclonal Antibody
Catalog #: AMRe08891




Background

Plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity. Plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.

Research Area

Image Data



Western blot analysis of Oligodendrocyte Specific Protein expression in SH-SY5Y cell lysate.

Note

For research use only.