

Product Name: Collagen VI (15E2) Rabbit Monoclonal Antibody**Catalog #: AMRe09217**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:100-1:200
Molecular Weight	147kDa

Antigen Information


Gene Name	
Alternative Names	Alpha 1 (VI) chain (61 AA); COL6A1; COL6A2; COL6A3; Collagen alpha 2(VI) chain; Collagen alpha 3(VI) chain; OPLL; PP3610;
Gene ID	
SwissProt ID	P12109/P12110/P12111
Immunogen	Recombinant protein of human Collagen VI

Background

Collagen VI acts as a cell-binding protein.

Research Area

Image Data



Western blot analysis of Collagen VI expression in Human brain lysate.