
Product Name: Cytokeratin 14 (12Q6) Rabbit Monoclonal Antibody**Catalog #: AMRe09724**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:2000-1:20000,IHC 1:200-1:500,ICC/IF 1:100-1:500,FC 1:50-1:100
Molecular Weight	52 kDa

Antigen Information


Gene Name	KRT14
Alternative Names	
Gene ID	3861.0
SwissProt ID	P02533
Immunogen	A synthetic peptide of human Cytokeratin 14

Background

Research Area

Signal Transduction

Image Data



Western blot analysis of Cytokeratin 14 in Mouse skin lysates using Cytokeratin 14 antibody.