

---

**Product Name: D4 GDI (13R12) Rabbit Monoclonal Antibody****Catalog #: AMRe09769**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,FC
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.5mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
<b>Purification</b>	Affinity purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:50-1:100,FC 1:50-1:100
<b>Molecular Weight</b>	23kDa

**Antigen Information**


<b>Gene Name</b>	ARHGDIB
<b>Alternative Names</b>	
<b>Gene ID</b>	397.0
<b>SwissProt ID</b>	P52566
<b>Immunogen</b>	A synthetic peptide of human D4 GDI

**Background**

## Research Area

Signal Transduction

## Image Data



Western blot detection of D4 GDI in Jurkat,U937 cell lysates using D4 GDI antibody(1:1000 diluted).