
Product Name: DAZL (11G2) Rabbit Monoclonal Antibody**Catalog #: AMRe09800**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.25mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:1000-1:5000,ICC/IF 1:200-1:500,IP 1:20-1:50
Molecular Weight	33kDa

Antigen Information


Gene Name	DAZL
Alternative Names	DAZ homolog; DAZ like autosomal; DAZH; DAZL; DAZL1; DAZLA; Deleted in azoospermia like 1; Germline specific RNA binding protein; SPGY like autosomal; SPGYLA; Tpx2;
Gene ID	1618.0
SwissProt ID	Q92904
Immunogen	Recombinant protein of human DAZL

Background

RNA-binding protein, which is essential for gametogenesis. Plays a central role during spermatogenesis. May act by binding to the 3'-UTR of mRNA and thereby regulating the translation of key transcripts. RNA-binding protein, which is essential for gametogenesis in both males and females. Plays a central role during spermatogenesis. Acts by binding to the 3'-UTR of mRNA, specifically recognizing GUU triplets, and thereby regulating the translation of key transcripts (By similarity).

Research Area

Image Data



Western blot analysis of DAZL expression in human testis lysate.