

Product Name: Dnmt1 (10B9) Rabbit Monoclonal Antibody**Catalog #: AMRe10087**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:200-1:1000,FC 1:50-1:200
Molecular Weight	183 kDa

Antigen Information


Gene Name	DNMT1
Alternative Names	
Gene ID	1786.0
SwissProt ID	P26358
Immunogen	A synthetic peptide of human Dnmt1

Background

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Dnmt1 in HEK293 , NIH/3T3, C6 lysates using Dnmt1 antibody.