
Product Name: DOPA Decarboxylase (18X12) Rabbit Monoclonal Antibody**Catalog #: AMRe10115**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	54kDa

Antigen Information


Gene Name	DDC
Alternative Names	AADC; Aromatic L Amino Acid Decarboxylase; DDC; DOPA decarboxylase (aromatic L-amino acid decarboxylase); DOPA decarboxylase;
Gene ID	1644.0
SwissProt ID	P20711
Immunogen	A synthetic peptide of human DOPA Decarboxylase

Background

Catalyzes the decarboxylation of L-3,4-dihydroxyphenylalanine (DOPA) to dopamine, L-5-hydroxytryptophan to serotonin and L-tryptophan to tryptamine. Catalyzes the decarboxylation of L-3,4-dihydroxyphenylalanine (DOPA) to dopamine, L-5-hydroxytryptophan to serotonin and L-tryptophan to tryptamine.

Research Area

Image Data



Western blot analysis of DOPA Decarboxylase expression in (1) 293T cell lysate; (2) RAW 264.7 cell lysate.