

Product Name: ERp57 (7K11) Rabbit Monoclonal Antibody**Catalog #: AMRe10613**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IF-P
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:50-1:100,IF-P 1:50-1:100
Molecular Weight	57kDa

Antigen Information

Gene Name	PDIA3
Alternative Names	ERp57; ERp60; ERp61; Glucose Regulated Protein 58 Kd; HsT17083; p58; PDIA3; Phospholipase C alpha; PI PLC;
Gene ID	2923.0
SwissProt ID	P30101
Immunogen	A synthetic peptide of human PDIA3


Background

Detected in the flagellum and head region of spermatozoa (at protein level).

Research Area

Antigen processing and presentation;

Image Data



Western blot analysis of ERp57 expression in HepG2 cell lysate.