

Product Name: ERp72 (5U2) Rabbit Monoclonal Antibody
Catalog #: AMRe10615

Summary

Production Name	ERp72 (5U2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Immunogen

Gene Name	PDIA4
Alternative Names	
Gene ID	9601.0
SwissProt ID	P13667. A synthetic peptide of human ERp72

Application


Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:50-1:100
Molecular Weight	73kDa

Background

Product Name: ERp72 (5U2) Rabbit Monoclonal Antibody
Catalog #: AMRe10615

Research Area

Image Data



Western blot detection of ERp72 in PC-3 cell lysates using ERp72 antibody(1:1000 diluted).

Note

For research use only.