

Product Name: ETS1 (11S7) Rabbit Monoclonal Antibody
Catalog #: AMRe10648

Summary

Production Name	ETS1 (11S7) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	ETS1
Alternative Names	C ets 1 protein; ETS 1; ETS1; ETS1 oncogene; EWSR2; Oncogene ETS1; P54; Protein C-ets-1;
Gene ID	2113.0
SwissProt ID	P14921.

Application

Dilution Ratio	WB 1:1000, IHC-P/IF-P 1:1000
Molecular Weight	50kDa

Product Name: ETS1 (11S7) Rabbit Monoclonal Antibody
Catalog #: AMRe10648




Background

Transcription factor. Transcription factor. Directly controls the expression of cytokine and chemokine genes in a wide variety of different cellular contexts. May control the differentiation, survival and proliferation of lymphoid cells. May also regulate angiogenesis through regulation of expression of genes controlling endothelial cell migration and invasion.

Research Area

Image Data



Western blot analysis of Insulin Receptor R expression in Daudi cell lysate.

Note

For research use only.