
Product Name: Fetuin A (6M11) Rabbit Monoclonal Antibody**Catalog #: AMRe10910**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	39kDa

Antigen Information

Gene Name	AHSG
Alternative Names	16G2; A2HS; AHS; Ahsg; FETUA; FETUB; FetuinA; FetuinB; Gugu; HSGA; IRL685; PRO2743;
Gene ID	197.0
SwissProt ID	P02765
Immunogen	A synthetic peptide of human Fetuin A

Background


Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and

barium ions. Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.

Research Area

Cardiovascular

Image Data



Western blot analysis of Fetuin A expression in Human plasma lysate.