

Product Name: Flotillin 2 (1C5) Rabbit Monoclonal Antibody
Catalog #: AMRe11028

Summary

Production Name	Flotillin 2 (1C5) Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	FLOT2
Alternative Names	epidermal surface antigen; ESA; flotillin-2; Membrane component chromosome 17 surface marker 1; ESA1; M17S1
Gene ID	2319.0
SwissProt ID	Q14254.

Application

Dilution Ratio	WB 1:500-1:2000,IP 1:20-1:50
Molecular Weight	47kDa

Product Name: Flotillin 2 (1C5) Rabbit Monoclonal Antibody
Catalog #: AMRe11028




Background

May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles. May be involved in epidermal cell adhesion and epidermal structure and function. May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles. May be involved in epidermal cell adhesion and epidermal structure and function.

Research Area

Insulin_Receptor;

Image Data



Western blot analysis of Flotillin-2 expression in (1) HeLa cell lysate;(2) 293 cell lysate.

Note

For research use only.