

Product Name: GCDFP 15 (8J14) Rabbit Monoclonal Antibody
Catalog #: AMRe11343

Summary

Production Name	GCDFP 15 (8J14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PIP
Alternative Names	BRST2; GCDGP15; gp17; GPIP4; PIP; Prolactin induced protein; SABP; Secretory actin binding protein;
Gene ID	5304.0
SwissProt ID	P12273.

Application

Dilution Ratio	WB 1:1000, IHC-P/IF-P 1:100-1:200, ICC/IF 1:50-1:100, FCM 1:20-1:50
Molecular Weight	17kDa

Product Name: GCDFP 15 (8J14) Rabbit Monoclonal Antibody
Catalog #: AMRe11343




Background

Belongs to the PIP family.

Research Area

Image Data



Western blot analysis of GCDFP 15 expression in T47D cell lysate.

Note

For research use only.