

Product Name: GM130 (6E7) Rabbit Monoclonal Antibody**Catalog #: AMRe11527**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.3mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:100-1:500,ICC/IF 1:50-1:200,FC 1:20-1:50,IP 1:20-1:50
Molecular Weight	113 kDa

Antigen Information

Gene Name	GOLGA2
Alternative Names	
Gene ID	2801.0
SwissProt ID	Q08379
Immunogen	A synthetic peptide of human GM130

Background

Research Area

Image Data

Western blot analysis of GM130 in HeLa, lysates using GM130 antibody.

