

Product Name: GRAF (2M19) Rabbit Monoclonal Antibody**Catalog #: AMRe11734**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:50-1:100
Molecular Weight	92kDa

Antigen Information


Gene Name	ARHGAP26
Alternative Names	GRAF; GRAF1; OPHN1L; OPHN1L1;
Gene ID	23092.0
SwissProt ID	Q9UNA1
Immunogen	A synthetic peptide of human GRAF

Background

GTPase-activating protein for RHOA and CDC42. GTPase-activating protein for RHOA and CDC42.

Research Area

Image Data



Western blot analysis of GRAF expression in HeLa cell lysate.