
Product Name: GSK3 alpha/beta (3Q8) Rabbit Monoclonal Antibody**Catalog #: AMRe11819**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:500-1:2000,ICC/IF 1:100-1:200,FC 1:50-1:100,IP 1:20-1:50
Molecular Weight	51,47kDa

Antigen Information


Gene Name	GSK3A
Alternative Names	
Gene ID	2931.0
SwissProt ID	P49840
Immunogen	Recombinant protein of human GSK3 alpha

Background

Research Area

Neuroscience

Image Data



Western blot analysis of extracts from A549 cells using GSK3 alpha/beta (3Q8) Rabbit Monoclonal Antibody at 1:1000.