

Product Name: HB9/HLXB9 (17W1) Rabbit Monoclonal Antibody
Catalog #: AMRe11908



Summary

Production Name	HB9/HLXB9 (17W1) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	MNX1
Alternative Names	HB9; HLXB9; Homeobox HB9; Homeobox protein HB9; HOXHB9; MNX1; Motor neuron and pancreas homeobox protein 1; SCRA1;
Gene ID	3110.0
SwissProt ID	P50219.

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	41kDa

Product Name: HB9/HLXB9 (17W1) Rabbit Monoclonal Antibody
Catalog #: AMRe11908




Background

Putative transcription factor involved in pancreas development and function. Putative transcription factor involved in pancreas development and function.

Research Area

Image Data



Western blot analysis of HB9/HLXB9 expression in Molt-4 cell lysate.

Note

For research use only.