

Product Name: HIRA (17J17) Rabbit Monoclonal Antibody
Catalog #: AMRe12038

Summary

Production Name	HIRA (17J17) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ICC/IF,FC,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	HIRA
Alternative Names	DGCR1; HIR; HIRA; TUP 1; TUPLE1;
Gene ID	7290.0
SwissProt ID	P54198.

Application

Dilution Ratio	WB 1:1000-1:5000, ICC/IF 1:50-1:200, FCM 1:100-1:1000, IP 1:10-1:100
Molecular Weight	112kDa


Product Name: HIRA (17J17) Rabbit Monoclonal Antibody
Catalog #: AMRe12038

Background

Cooperates with ASF1A to promote replication-independent chromatin assembly. Required for the periodic repression of histone gene transcription during the cell cycle. Cooperates with ASF1A to promote replication-independent chromatin assembly. Required for the periodic repression of histone gene transcription during the cell cycle. Required for the formation of senescence-associated heterochromatin foci (SAHF) and efficient senescence-associated cell cycle exit.

Research Area

Image Data



Western blot analysis of HIRA expression in HeLa cell lysate.

Note

For research use only.