

Product Name: Huntingtin (2P4) Rabbit Monoclonal Antibody
Catalog #: AMRe12284



Summary

Production Name	Huntingtin (2P4) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	HTT
Alternative Names	HD; HDH; HTT; Huntingtin; IT15; ZHD;
Gene ID	3064.0
SwissProt ID	P42858.

Application

Dilution Ratio	WB 1:1000, IHC-P/IF-P 1:100, ICC/IF 1:200, FCM 1:50
Molecular Weight	348kDa

Product Name: Huntingtin (2P4) Rabbit Monoclonal Antibody
Catalog #: AMRe12284




Background

May play a role in microtubule-mediated transport or vesicle function. [Huntingtin]: May play a role in microtubule-mediated transport or vesicle function.

Research Area

Image Data



Western blot analysis of Huntingtin expression in HeLa cell lysate.

Note

For research use only.