

Product Name: ICAD (16K15) Rabbit Monoclonal Antibody**Catalog #: AMRe12333**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.33mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:200-1:1000,ICC/IF 1:200-1:500
Molecular Weight	37kDa

Antigen Information

Gene Name	DFFA
Alternative Names	DFF1; DFF35; DFF45; DFFA; H13; ICAD; ICAD L; ICAD S;
Gene ID	1676.0
SwissProt ID	O00273
Immunogen	A synthetic peptide of human ICAD


Background

Inhibitor of the caspase-activated DNase (DFF40). Inhibitor of the caspase-activated DNase (DFF40).

Research Area

Cell Biology

Image Data



Western blot analysis of ICAD expression in HeLa cell lysate.