
Product Name: IKK alpha/beta (10K2) Rabbit Monoclonal Antibody**Catalog #: AMRe12472**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IP 1:20-1:50
Molecular Weight	85 kDa

Antigen Information


Gene Name	CHUK/IKBKB
Alternative Names	
Gene ID	1147/3551
SwissProt ID	O15111/O14920
Immunogen	Recombinant protein of human IKK alpha/beta

Background

Research Area

Signal Transduction

Image Data



Western blot analysis of IKK alpha/beta in 293T lysates using IKK alpha/beta antibody.