
Product Name: IL18 binding protein (18K17) Rabbit Monoclonal Antibody**Catalog #: AMRe12513**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	21kDa

Antigen Information


Gene Name	IL18BP
Alternative Names	IL18 binding protein; IL18 Bpa; IL18BP; IL18BPa; Interleukin 18 binding protein; MC51L 53L 54L homolog gene product; Tadekinig alfa;
Gene ID	10068.0
SwissProt ID	O95998
Immunogen	A synthetic peptide of human IL18 binding protein

Background

Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response. Isoform A binds to IL-18 and inhibits its activity. Functions as an inhibitor of the early TH1 cytokine response.

Research Area

Image Data



Western blot analysis of IL18 binding protein in HeLa cell lysate.