

Product Name: IL9 (14B18) Rabbit Monoclonal Antibody**Catalog #: AMRe12574**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, FC, IP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000, ICC/IF 1:200-1:500, FC 1:200-1:500, IP 1:20-1:50
Molecular Weight	16kDa

Antigen Information

Gene Name	IL9
Alternative Names	HP40; IL9; Interleukin 9 MCGF; P40; TCGF 3;
Gene ID	3578.0
SwissProt ID	P15248
Immunogen	Recombinant protein of human IL-9


Background

Supports IL-2 independent and IL-4 independent growth of helper T-cells. Supports IL-2 independent and IL-4 independent

growth of helper T-cells.

Research Area

Image Data



Western blot analysis of IL9 in Recombinant Human IL9 protein cell lysate.