
Product Name: Intra Acrosomal Protein (7H17) Rabbit Monoclonal Antibody**Catalog #: AMRe12688**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	28kDa

Antigen Information


Gene Name	ACRV1
Alternative Names	ACRV1; Msa63; Sp10; SPACA2;
Gene ID	56.0
SwissProt ID	P26436
Immunogen	A synthetic peptide of human Intra Acrosomal Protein

Background

N/A

Research Area

Image Data



Western blot analysis of Intra Acrosomal Protein expression in Human testis lysate.