

Product Name: IRS2 (17U8) Rabbit Monoclonal Antibody
Catalog #: AMRe12760

Summary

Production Name	IRS2 (17U8) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	IRS2
Alternative Names	Insulin receptor substrate 2; IRS 2; IRS-2;
Gene ID	8660.0
SwissProt ID	Q9Y4H2.

Application

Dilution Ratio	WB 1:1000-1:5000
Molecular Weight	137kDa


Background

Product Name: IRS2 (17U8) Rabbit Monoclonal Antibody
Catalog #: AMRe12760

May mediate the control of various cellular processes by insulin. May mediate the control of various cellular processes by insulin.

Research Area

Image Data



Western blot analysis of IRS2 expression in HEK293 cell treated with insulin.

Note

For research use only.