
Product Name: JNK (19C18) Rabbit Monoclonal Antibody**Catalog #: AMRe12841**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF,FC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ICC/IF 1:500-1:1000,FC 1:50-1:200,IP 1:20-1:50
Molecular Weight	48,53kDa

Antigen Information


Gene Name	MAPK8
Alternative Names	
Gene ID	5599.0
SwissProt ID	P45983
Immunogen	Recombinant protein of human JNK1

Background

Research Area

Signal Transduction

Image Data



Western blot detection of JNK1/2/3 in Jurkat,Rat Brain,C6,CHO-K1,HeLa cell lysates using JNK1/2/3 antibody(1:1000 diluted).