

**Product Name: KLK2 (11A15) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe13075**

---

## Summary

|                        |                                         |
|------------------------|-----------------------------------------|
| <b>Production Name</b> | KLK2 (11A15) Rabbit Monoclonal Antibody |
| <b>Description</b>     | Recombinant rabbit monoclonal antibody  |
| <b>Host</b>            | Rabbit                                  |
| <b>Application</b>     | WB, ICC/IF                              |
| <b>Reactivity</b>      | Human, Mouse, Rat                       |

## Performance

|                     |                                                                                                                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                                                                                                                 |
| <b>Modification</b> | Unmodified                                                                                                                                                                                   |
| <b>Isotype</b>      | IgG                                                                                                                                                                                          |
| <b>Clonality</b>    | Monoclonal                                                                                                                                                                                   |
| <b>Form</b>         | Liquid                                                                                                                                                                                       |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                                                     |
| <b>Buffer</b>       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |
| <b>Purification</b> | Affinity purification                                                                                                                                                                        |

## Immunogen

|                          |                     |
|--------------------------|---------------------|
| <b>Gene Name</b>         | KLK2                |
| <b>Alternative Names</b> | hK2; hGK-1; KLK2A2; |
| <b>Gene ID</b>           | 3817.0              |
| <b>SwissProt ID</b>      | P20151.             |

## Application

|                         |                                     |
|-------------------------|-------------------------------------|
| <b>Dilution Ratio</b>   | WB 1:1000-1:5000, ICC/IF 1:50-1:100 |
| <b>Molecular Weight</b> | 29kDa                               |

**Product Name: KLK2 (11A15) Rabbit Monoclonal Antibody**  
**Catalog #: AMRe13075**


---

## Background

Glandular kallikreins cleave Met-Lys and Arg-Ser bonds in kininogen to release Lys-bradykinin. Glandular kallikreins cleave Met-Lys and Arg-Ser bonds in kininogen to release Lys-bradykinin.

## Research Area

## Image Data



Western blot analysis of KLK2 expression in MOLT4 cell lysate.

## Note

For research use only.