

Product Name: Laminin gamma 1 (9A17) Rabbit Monoclonal Antibody
Catalog #: AMRe13195

Summary

Production Name	Laminin gamma 1 (9A17) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,FC,IP,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Immunogen

Gene Name	LAMB2; lamc1; Laminin B2 chain; Laminin subunit gamma-1; Laminin-1 subunit gamma; Laminin-10 subunit gamma; Laminin-11 subunit gamma; Laminin-2 subunit gamma; Laminin-3 subunit gamma; Laminin-4 subunit gamma; Laminin-6 subunit gamma; Laminin-7 subunit gamma; Laminin-8 subunit gamma; Laminin-9 subunit gamma;
Alternative Names	
Gene ID	
SwissProt ID	P11047/Q5VYE7/Q6NVY8. Recombinant protein of human Laminin gamma 1

Application

Dilution Ratio	WB 1:1000, IHC-P/IF-P 1:500, FCM 1:20, IP 1:20
-----------------------	--

**Product Name: Laminin gamma 1 (9A17) Rabbit
Monoclonal Antibody
Catalog #: AMRe13195**




Molecular Weight 178kDa

Background

Binding to cells via a high affinity receptor, laminin is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components.

Research Area

Image Data



Western blot detection of Laminin gamma 1 in HeLa, MCF-7, A549, HCT116 cell lysates using Laminin gamma 1 antibody (1:1000 diluted).

Note

For research use only.