

Product Name: LDHB (13O6) Rabbit Monoclonal Antibody
Catalog #: AMRe13264

Summary

Production Name	LDHB (13O6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	LDHB
Alternative Names	HEL S 281; LDH B; LDH H; LDH-B; LDH-H; LDHB; LDHBD; LDHH; TRG5;
Gene ID	3945.0
SwissProt ID	P07195.

Application

Dilution Ratio	WB 1:5000-1:50000, IHC-P/IF-P 1:100-1:200
Molecular Weight	37kDa

Product Name: LDHB (13O6) Rabbit Monoclonal Antibody
Catalog #: AMRe13264




Background

Fermentation; pyruvate fermentation to lactate; (S)-lactate from pyruvate: step 1/1.

Research Area

Image Data



Western blot analysis of LDHB expression in HeLa cell lysate.

Note

For research use only.