

Product Name: LGALS3BP (7X1) Rabbit Monoclonal Antibody**Catalog #: AMRe13287**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,FC 1:50-1:100
Molecular Weight	65kDa

Antigen Information


Gene Name	LGALS3BP
Alternative Names	
Gene ID	3959.0
SwissProt ID	Q08380
Immunogen	Recombinant protein of human LGALS3BP

Background

Research Area

Cancer

Image Data



Western blot detection of LGALS3BP in HT1080, HepG2, A549 cell lysates using LGALS3BP antibody(1:1000 diluted).