

Product Name: LMO4 (7B8) Rabbit Monoclonal Antibody
Catalog #: AMRe13358

Summary

Production Name	LMO4 (7B8) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	LMO4
Alternative Names	LIM domain only protein 4; LMO4;;
Gene ID	8543.0
SwissProt ID	P61968.

Application

Dilution Ratio	WB 1:1000-1:5000, IP 1:20-1:100
Molecular Weight	18kDa


Product Name: LMO4 (7B8) Rabbit Monoclonal Antibody
Catalog #: AMRe13358

Background

Probable transcriptional factor. Probable transcriptional factor.

Research Area

Image Data



Western blot analysis of LMO4 expression in Jurkat cell lysate.

Note

For research use only.