
Product Name: LPCAT1 (13M5) Rabbit Monoclonal Antibody**Catalog #: AMRe13387**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:2000,ICC/IF 1:100-1:200,FC 1:20-1:50,IP 1:20-1:50
Molecular Weight	59kDa


Antigen Information

Gene Name	LPCAT1
Alternative Names	
Gene ID	79888.0
SwissProt ID	Q8NF37
Immunogen	Recombinant protein of human LPCAT1

Background

Research Area

Image Data



Western blot detection of LPCAT1 in A549 cell lysates using LPCAT1 antibody(1:1000 diluted).