

Product Name: LRG1 (8F15) Rabbit Monoclonal Antibody**Catalog #: AMRe13418**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP,IF-P
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:100,IP 1:20-1:50,IF-P 1:50-1:100
Molecular Weight	38kDa

Antigen Information

Gene Name	LRG1
Alternative Names	A2GL; HMFT1766; LRG; LRG1;
Gene ID	116844.0
SwissProt ID	P02750
Immunogen	A synthetic peptide of human LRG1


Background

Transforming growth factor beta receptor binding.

Research Area

Immunology

Image Data



Western blot analysis of LRG1 expression in Human plasma lysate.