

Product Name: Mammaglobin A (8R14) Rabbit Monoclonal Antibody**Catalog #: AMRe13613**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IF-P
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:200-1:1000,IF-P 1:200-1:1000
Molecular Weight	10kDa

Antigen Information


Gene Name	SCGB2A2
Alternative Names	Mammaglobin 1; Mammaglobin-A; MGB1; Scgb2a2; SG2A2; UGB2;
Gene ID	4250.0
SwissProt ID	Q13296
Immunogen	Recombinant protein of human Mammaglobin A

Background

Belongs to the secretoglobin family. Lipophilin subfamily.

Research Area

Image Data



Western blot analysis of Mammaglobin A expression in Human breast cancer lysate.