

**Product Name: MTX3 (10I6) Rabbit Monoclonal Antibody****Catalog #: AMRe14224**

For research use only.

**Summary**

|                      |                                                                                                                             |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>   | Recombinant rabbit monoclonal antibody                                                                                      |
| <b>Host</b>          | Rabbit                                                                                                                      |
| <b>Application</b>   | WB                                                                                                                          |
| <b>Reactivity</b>    | Human, Mouse, Rat                                                                                                           |
| <b>Conjugation</b>   | Unconjugated                                                                                                                |
| <b>Modification</b>  | Unmodified                                                                                                                  |
| <b>Isotype</b>       | IgG                                                                                                                         |
| <b>Clonality</b>     | Monoclonal                                                                                                                  |
| <b>Form</b>          | Liquid                                                                                                                      |
| <b>Concentration</b> | 0.5mg/ml. The concentration of this product may be batch-dependent.                                                         |
| <b>Storage</b>       | Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.                                                 |
| <b>Shipping</b>      | Ice bags                                                                                                                    |
| <b>Buffer</b>        | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein. |
| <b>Purification</b>  | Affinity purification                                                                                                       |

**Application**

|                         |                  |
|-------------------------|------------------|
| <b>Dilution Ratio</b>   | WB 1:1000-1:5000 |
| <b>Molecular Weight</b> | 35kDa            |

**Antigen Information**


|                          |                                   |
|--------------------------|-----------------------------------|
| <b>Gene Name</b>         | MTX3                              |
| <b>Alternative Names</b> |                                   |
| <b>Gene ID</b>           | 345778.0                          |
| <b>SwissProt ID</b>      | Q5HYI7                            |
| <b>Immunogen</b>         | A synthetic peptide of human MTX3 |

**Background**

## Research Area

Signal Transduction

## Image Data



Western blot detection of MTX3 in mouse brain cell lysates using MTX3 antibody(1:1000 diluted).