
Product Name: MYH11 (10S12) Rabbit Monoclonal Antibody**Catalog #: AMRe14289**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:200-1:1000,ICC/IF 1:20-1:50,FC 1:20-1:50
Molecular Weight	227kDa

Antigen Information


Gene Name	MYH11
Alternative Names	AAT4; FAA4; MYH11; SM1; SM2; SMHC; SMMHC;
Gene ID	4629.0
SwissProt ID	P35749
Immunogen	A synthetic peptide of human smooth muscle Myosin heavy chain 11

Background

Muscle contraction. Muscle contraction.

Research Area

Image Data



Western blot analysis of Myosin, smooth muscle heavy chain 1 and 2 expression in Human bladder lysate.