

Product Name: MYH11 (6C13) Rabbit Monoclonal Antibody
Catalog #: AMRe14290

Summary

Production Name	MYH11 (6C13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Immunogen

Gene Name	MYH11
Alternative Names	AAT4; FAA4; MYH11; Myosin 11; Myosin heavy chain 11; Myosin heavy chain smooth muscle isoform; Myosin heavy polypeptide 11 smooth muscle; SMHC; SMMHC;
Gene ID	4629.0
SwissProt ID	P35749. A synthetic peptide of human Myosin, smooth muscle heavy chain 1 and 2

Application

Dilution Ratio	WB 1:2000-1:20000, IHC-P/IF-P 1:100-1:500
Molecular Weight	227kDa

Product Name: MYH11 (6C13) Rabbit Monoclonal Antibody
Catalog #: AMRe14290




Background

Muscle contraction. Muscle contraction.

Research Area

Image Data



Western blot analysis of extracts from Mouse spleen tissue using RM5298 at 1:1000.

Note

For research use only.