

Product Name: NeuN (11E16) Rabbit Monoclonal Antibody**Catalog #: AMRe14584**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.43mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:500-1:2000,ICC/IF 1:100-1:200,FC 1:50-1:200
Molecular Weight	34kDa

Antigen Information


Gene Name	RBFOX3
Alternative Names	
Gene ID	146713.0
SwissProt ID	A6NFN3
Immunogen	A synthetic peptide of human NeuN

Background

Research Area

Neuroscience

Image Data



Western blot detection of NeuN in Rat Brain lysates using NeuN antibody(1:1000 diluted).