

Product Name: Oct-2 (16S17) Rabbit Monoclonal Antibody**Catalog #: AMRe15101**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IF-P
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:100-1:200,IF-P 1:100-1:200
Molecular Weight	51kDa

Antigen Information


Gene Name	POU2F2
Alternative Names	
Gene ID	5452.0
SwissProt ID	P09086
Immunogen	Recombinant protein of human Oct-2

Background

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of extracts from Jurkat cells using Oct-2 (16S17) Rabbit Monoclonal Antibody at 1:1000.