

Product Name: P Glycoprotein (3J14) Rabbit Monoclonal Antibody
Catalog #: AMRe15561

Summary

Production Name	P Glycoprotein (3J14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
Purification	Affinity purification

Immunogen

Gene Name	ABCB1
Alternative Names	
Gene ID	5243.0
SwissProt ID	P08183.

Application


Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:200-1:1000
Molecular Weight	141 kDa

Background

Product Name: P Glycoprotein (3J14) Rabbit Monoclonal Antibody
Catalog #: AMRe15561

Research Area

Image Data



Western blot analysis of P Glycoprotein in HepG2, RAW 264.7, C6 lysates using P Glycoprotein antibody.

Note

For research use only.