

Product Name: P4HB (15Y14) Rabbit Monoclonal Antibody**Catalog #: AMRe15635**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.37mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:200-1:1000,ICC/IF 1:100-1:500,FC 1:100-1:1000
Molecular Weight	57kDa

Antigen Information


Gene Name	P4HB
Alternative Names	
Gene ID	5034.0
SwissProt ID	P07237
Immunogen	A synthetic peptide of human P4HB

Background

Research Area

Tags & Cell Markers

Image Data



Western blot detection of P4HB in HCT116, Hela, HepG2, Jurkat, MCF-7 cell lysates using P4HB antibody (1:1000 diluted).