

Product Name: Parkin (17H13) Rabbit Monoclonal Antibody
Catalog #: AMRe15758

Summary

Production Name	Parkin (17H13) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PRKN
Alternative Names	AR-JP; LPRS2; PARK2; parkin; parkin 2; PDJ; PRKN; PRKN2;
Gene ID	5071.0
SwissProt ID	O60260.

Application

Dilution Ratio	WB 1:1000
Molecular Weight	52kDa

Product Name: Parkin (17H13) Rabbit Monoclonal Antibody
Catalog #: AMRe15758



Background

Parkin is a zinc-finger protein that is related to ubiquitin at the amino terminus. May protect neurons against alpha synuclein toxicity, proteasomal dysfunction, GPR37 accumulation, and kainate-induced excitotoxicity. May play a role in controlling neurotransmitter trafficking at the presynaptic terminal and in calcium-dependent exocytosis. Regulates cyclin-E during neuronal apoptosis. Functions within a multiprotein E3 ubiquitin ligase complex, catalyzing the covalent attachment of ubiquitin moieties onto substrate proteins (PubMed:< a href="http://www.uniprot.org/citations/10888878"

target="_blank">10888878, PubMed:< a href="http://www.uniprot.org/citations/10973942"

target="_blank">10973942, PubMed:< a href="http://www.uniprot.org/citations/11431533"

target="_blank">11431533, PubMed:< a href="http://www.uniprot.org/citations/12150907"

target="_blank">12150907, PubMed:< a href="http://www.uniprot.org/citations/12628165"

target="_blank">12628165, PubMed:< a href="http://www.uniprot.org/citations/15105460"

target="_blank">15105460, PubMed:< a href="http://www.uniprot.org/citations/16135753"

target="_blank">16135753, PubMed:< a href="http://www.uniprot.org/citations/21376232"

target="_blank">21376232, PubMed:< a href="http://www.uniprot.org/citations/21532592"

target="_blank">21532592, PubMed:< a href="http://www.uniprot.org/citations/23754282"

target="_blank">23754282, PubMed:< a href="http://www.uniprot.org/citations/23620051"

target="_blank">23620051, PubMed:< a href="http://www.uniprot.org/citations/24660806"

target="_blank">24660806, PubMed:< a href="http://www.uniprot.org/citations/24751536"

target="_blank">24751536, PubMed:< a href="http://www.uniprot.org/citations/32047033"

target="_blank">32047033, PubMed:< a href="http://www.uniprot.org/citations/29311685"

target="_blank">29311685, PubMed:< a href="http://www.uniprot.org/citations/22396657"

target="_blank">22396657). Substrates include SYT11 and VDAC1 (PubMed:< a

href="http://www.uniprot.org/citations/32047033" target="_blank">32047033, PubMed:< a

href="http://www.uniprot.org/citations/29311685" target="_blank">29311685). Other substrates are BCL2, CCNE1,

GPR37, RHOT1/MIRO1, MFN1, MFN2, STUB1, SNCAIP, SEPTIN5, TOMM20, USP30, ZNF746, MIRO1 and AIMP2 (PubMed:< a

href="http://www.uniprot.org/citations/10888878" target="_blank">10888878, PubMed:< a

href="http://www.uniprot.org/citations/10973942" target="_blank">10973942, PubMed:< a

href="http://www.uniprot.org/citations/11431533" target="_blank">11431533, PubMed:< a

href="http://www.uniprot.org/citations/12150907" target="_blank">12150907, PubMed:< a

href="http://www.uniprot.org/citations/12628165" target="_blank">12628165, PubMed:< a

href="http://www.uniprot.org/citations/15105460" target="_blank">15105460, PubMed:< a

href="http://www.uniprot.org/citations/16135753" target="_blank">16135753, PubMed:< a

href="http://www.uniprot.org/citations/21376232" target="_blank">21376232, PubMed:< a

href="http://www.uniprot.org/citations/21532592" target="_blank">21532592, PubMed:< a

href="http://www.uniprot.org/citations/23754282" target="_blank">23754282, PubMed:< a

Product Name: Parkin (17H13) Rabbit Monoclonal Antibody
Catalog #: AMRe15758



<http://www.uniprot.org/citations/23620051> target="_blank">23620051, PubMed:24660806, PubMed:24751536, PubMed:22396657). Mediates monoubiquitination as well as 'Lys-6', 'Lys-11', 'Lys-48'-linked and 'Lys-63'-linked polyubiquitination of substrates depending on the context (PubMed:19229105, PubMed:20889974, PubMed:25621951, PubMed:32047033, PubMed:25474007). Participates in the removal and/or detoxification of abnormally folded or damaged protein by mediating 'Lys-63'-linked polyubiquitination of misfolded proteins such as PARK7: 'Lys-63'-linked polyubiquitinated misfolded proteins are then recognized by HDAC6, leading to their recruitment to aggresomes, followed by degradation (PubMed:17846173, PubMed:19229105). Mediates 'Lys-63'-linked polyubiquitination of a 22 kDa O-linked glycosylated isoform of SNCAIP, possibly playing a role in Lewy-body formation (PubMed:11431533, PubMed:11590439, PubMed:15105460, PubMed:19229105, PubMed:15728840). Mediates monoubiquitination of BCL2, thereby acting as a positive regulator of autophagy (PubMed:20889974). Protects against mitochondrial dysfunction during cellular stress, by acting downstream of PINK1 to coordinate mitochondrial quality control mechanisms that remove and replace dysfunctional mitochondrial components (PubMed:32047033, PubMed:19029340, PubMed:19966284, PubMed:23620051, PubMed:24896179, PubMed:25527291, PubMed:18957282, PubMed:21376232, PubMed:22396657, PubMed:24660806, PubMed:25474007, PubMed:<a

Web: <https://www.enkilife.com> E-mail: order@enkilife.com techsupport@enkilife.com Tel: 0086-27-87002838


Product Name: Parkin (17H13) Rabbit Monoclonal Antibody
Catalog #: AMRe15758

[25527291](http://www.uniprot.org/citations/25527291) , PubMed:22396657 , PubMed:23933751). Preferentially assembles 'Lys-6', 'Lys-11'- and 'Lys-63'-linked polyubiquitin chains, leading to mitophagy (PubMed:25621951 , PubMed:32047033). The PINK1-PRKN pathway also promotes fission of damaged mitochondria by PINK1-mediated phosphorylation which promotes the PRKN-dependent degradation of mitochondrial proteins involved in fission such as MFN2 (PubMed:23620051). This prevents the refusion of unhealthy mitochondria with the mitochondrial network or initiates mitochondrial fragmentation facilitating their later engulfment by autophagosomes (PubMed:23620051). Regulates motility of damaged mitochondria via the ubiquitination and subsequent degradation of MIRO1 and MIRO2; in motor neurons, this likely inhibits mitochondrial intracellular anterograde transport along the axons which probably increases the chance of the mitochondria undergoing mitophagy in the soma (PubMed:22396657). Involved in mitochondrial biogenesis via the 'Lys-48'-linked polyubiquitination of transcriptional repressor ZNF746/PARIS which leads to its subsequent proteasomal degradation and allows activation of the transcription factor PPARGC1A (PubMed:21376232). Limits the production of reactive oxygen species (ROS) (PubMed:18541373). Regulates cyclin-E during neuronal apoptosis (PubMed:12628165). In collaboration with CHPF isoform 2, may enhance cell viability and protect cells from oxidative stress (PubMed:22082830). Independently of its ubiquitin ligase activity, protects from apoptosis by the transcriptional repression of p53/TP53 (PubMed:19801972). May protect neurons against alpha synuclein toxicity, proteasomal dysfunction, GPR37 accumulation, and kainate-induced excitotoxicity (PubMed:11439185). May play a role in controlling neurotransmitter trafficking at the presynaptic terminal and in calcium-dependent exocytosis. May represent a tumor suppressor gene (PubMed:12719539).

Research Area

Image Data

Product Name: Parkin (17H13) Rabbit Monoclonal Antibody
Catalog #: AMRe15758



Western blot analysis of Parkin expression in Jurkat cell lysate.

Note

For research use only.