

Product Name: PHLDA1 (11J6) Rabbit Monoclonal Antibody**Catalog #: AMRe16077**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:2000-1:20000,IHC 1:50-1:100,ICC/IF 1:100-1:200
Molecular Weight	45kDa

Antigen Information


Gene Name	PHLDA1
Alternative Names	
Gene ID	22822.0
SwissProt ID	Q8WV24
Immunogen	A synthetic peptide of human PHLDA1

Background

Research Area

Cell Biology

Image Data



Western blot detection of PHLDA1 in A375 cell lysates using PHLDA1 antibody(1:1000 diluted).