

**Product Name: Prolactin Receptor (2K9) Rabbit
Monoclonal Antibody
Catalog #: AMRe16522**

Summary

Production Name	Prolactin Receptor (2K9) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,IF-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type
Buffer	preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	PRLR
Alternative Names	CPRLP; HPRL; hPRL receptor; Lactogen receptor; MFAB; OPR; Pr-1; Pr-3; PRLR; Prolactin receptor a;
Gene ID	5618.0
SwissProt ID	P16471.

Application

Dilution Ratio	WB 1:1000-1:5000, IHC-P/IF-P 1:50-1:100, ICC/IF 1:100-1:200
Molecular Weight	70kDa

**Product Name: Prolactin Receptor (2K9) Rabbit
Monoclonal Antibody
Catalog #: AMRe16522**




Background

Activate multiple signal transduction pathways in breast cancer cells, including the Jak/Stat, PI3K/Akt, and MAPK pathways, leading to both pro-proliferative and anti-apoptotic downstream effects. This is a receptor for the anterior pituitary hormone prolactin (PRL). Acts as a prosurvival factor for spermatozoa by inhibiting sperm capacitation through suppression of SRC kinase activation and stimulation of AKT. Isoform 4 is unable to transduce prolactin signaling. Isoform 6 is unable to transduce prolactin signaling.

Research Area

Image Data



Western blot analysis of Prolactin Receptor expression in Jurkat cell lysate.

Note

For research use only.