

Product Name: Prostein (16Q10) Rabbit Monoclonal Antibody
Catalog #: AMRe16529

Summary

Production Name	Prostein (16Q10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	IHC-P,IF-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	SLC45A3
Alternative Names	IPCA6; PCANAP2; PCANAP6; PCANAP8; Prostein; PRST; Slc45a3;
Gene ID	85414.0
SwissProt ID	Q96JT2.

Application

Dilution Ratio	IHC-P/IF-P 1:200-1:2000
Molecular Weight	59kDa

Product Name: Prostein (16Q10) Rabbit Monoclonal Antibody
Catalog #: AMRe16529




Background

Marker for prostate cells. May be used, in case of prostate cancers, as a target antigen for prostate carcinomas-directed cytotoxic T-cell lymphocytes.

Research Area

Image Data



Western blot analysis of Prostein expression in LnCaP cell lysate.

Note

For research use only.