

Product Name: Pulmonary Surfactant Associated Protein D (13G17) Rabbit Monoclonal Antibody

Catalog #: AMRe16689

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IF-P
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:1000-1:5000,IF-P 1:1000-1:5000
Molecular Weight	38kDa

Antigen Information


Gene Name	Sftpd
Alternative Names	
Gene ID	20390.0
SwissProt ID	P50404
Immunogen	Recombinant protein of mouse Surfactant protein D

Background

Research Area

Immunology

Image Data



Western blot detection of Surfactant protein D in Rat lung cell lysates using Surfactant protein D antibody(1:1000 diluted).