

Product Name: RPS6 (6I7) Rabbit Monoclonal Antibody
Catalog #: AMRe17373



Summary

Production Name	RPS6 (6I7) Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
Purification	Affinity purification

Immunogen

Gene Name	RPS6
Alternative Names	
Gene ID	6194.0
SwissProt ID	P62753.A synthetic peptide of human RPS6

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:1000,ICC/IF 1:200-1:500,FC 1:200-1:500
Molecular Weight	29kDa

Background


Product Name: RPS6 (6I7) Rabbit Monoclonal Antibody
Catalog #: AMRe17373



Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot detection of RPS6 in K562,C6,3T3 cell lysates using RPS6 antibody(1:1000 diluted).

Note

For research use only.