

Product Name: SAA4 (11R1) Rabbit Monoclonal Antibody**Catalog #: AMRe17567**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:100,FC 1:20-1:50
Molecular Weight	15kDa

Antigen Information


Gene Name	SAA4
Alternative Names	CSAA; Saa4;
Gene ID	6291.0
SwissProt ID	P35542
Immunogen	A synthetic peptide of human SAA4

Background

Major acute phase reactant. Apolipoprotein of the HDL complex. Major acute phase reactant.

Research Area

Image Data



Western blot analysis of SAA4 expression in Human plasma cell lysate.